

Liberty

Dual Mobility Cup System



Liberty

Dual Mobility Cup System

Product Specifications



LIBERTY SC



LIBERTY AC



LIBERTY T

FIXATION	Cementless 0.8mm equatorial press-fit 5 fins	Cemented Grooves to facilitate cement fixation	Cementless Addition of flange with screwhole and two spikes
SIZES	44-68 mm		
SHELL MATERIAL	Stainless Steel (ISO 5832-9)		
LINER OPTIONS	ID 22mm (OD 44-46) and ID 28mm (OD 48-68)		
LINER MATERIAL	Ultra High Molecular Weight Polyethylene (UHMWPE)		
SHELL COATING	Titanium (300 µm) and HA Coating (100 µm)	N/A	Titanium (300 µm) and HA Coating (100 µm)

- The LIBERTY Dual Mobility Acetabular System design is based on the traditional dual articulation concept, which has been in clinical use since 1977.
- The LIBERTY SC and AC Dual Mobility Acetabular Systems were marketed in 2005 and have not been modified since then.
- The LIBERTY T option was first placed on the market in 2007 and has not been modified since then.
- Screws are available in 5mm diameter, length from 30 to 60mm in 5mm increments.



MicroPort Orthopedics Inc.
5677 Airline Road
Arlington, TN USA 38002
866 872 0211
ortho.microport.com



MicroPort Orthopedics BV
Hoogoorddreef 5
1101 BA Amsterdam
The Netherlands
+31 20 545 01 00

The CE-Marking of Conformity is applied per catalog number and appears on the outer package label, if applicable.
Trademarks and Registered marks of MicroPort Orthopedics Inc. © 2015 MicroPort Orthopedics Inc. Liberty is a trademark of ATR. All Rights Reserved. 010781